

Cross-Reactivity for Beckman AU Urine Drug Screen Analytes

Amphetamine

Compounds	Concentration equivalent to 1000 ng/mL Cutoff	% Cross-Reactivity
d-Amphetamine	1000	100
Benzphetamine*	1000	100
d,l-Methamphetamine	2100	47.6
d,l-Amphetamine	2150	46.5
l-Methamphetamine	3650	27.4
MDA	6500	15.4
Donepezil	11200	8.9
l-Amphetamine	11500	8.7
4-Chloramphetamine	12200	8.2
Phenmetrazine	13000	7.7
1,3 Dimethylpentylamine	14900	6.7
d,l-4-Methylamphetamine	16500	6.1
Phentermine	25000	4.0
MDEA	27200	3.7
MDMA	34300	2.9
Isometheptene	56000	1.8
Mephentermine	60000	1.7
erythro-Dihydrobupropion	82000	1.2
Fenfluramine	150000	0.7
Nor-pseudoephedrine	170000	0.6
Tranlycypromine	200000	0.5
Methoxyphenamine	360000	0.3
Propranolol	500000	0.2
Tyramine	600000	0.2
Phenylpropanolamine	2000000	0.05
l-Ephedrine	3500000	0.03
Bupropion	4500000	0.02
Chloroquine	4500000	0.02
d,l-Pseudoephedrine	8300000	0.01
Quinacrine	16500000	<0.01
Benzphetamine metabolizes to amphetamine and methamphetamine		
Selegiline metabolizes to l-amphetamine and l-methamphetamine		

Barbiturates

Compounds	Concentration equivalent to 200 ng/mL Cutoff	% Cross-Reactivity
Talbital	194	103.1
Secobarbital	200	100.0
Pentobarbital	252	79.4
Butabarbital	274	73.0
Aprobarbital	275	72.7
Alphenal	284	70.4
Butalbital	304	65.8
Cyclopentobarbital	304	65.8
Amobarbital	348	57.5
Butobarbital	349	57.3
Phenobarbital	509	39.3
5-Ethyl-5-(4-hydroxyphenyl) barbituric acid	927	21.6
Barbital	1278	15.6
Thiopental	28200	0.7

Buprenorphine

Compounds	Concentration equivalent to 5 ng/mL Cutoff	% Cross-Reactivity
Norbuprenorphine	4.6	92.6
Buprenorphine Glucuronide	1000	0.1
Norbuprenorphine Glucuronide	1000	0.1

Benzodiazepine

Compounds	Concentration equivalent to 200 ng/mL Cutoff	% Cross-Reactivity
Alprazolam	65	307.7
Diazepam	70	285.7
Tetrazepam	70	285.7
Nitrazepam	78	256.4
Estazolam	90	222.2
Phenazepam	90	222.2
Prazepam	90	222.2
α -Hydroxyalprazolam	100	200.0
Ketazolam	100	200.0
N-Desmethyldiazepam	110	181.8
Halazepam	110	181.8
α -Hydroxyalprazolam glucuronide†	110	181.8
N-Desalkylflurazepam	130	153.8
α -Hydroxytriazolam	130	153.8
Midazolam	130	153.8
Triazolam	130	153.8
Flunitrazepam	140	142.9
Temazepam	140	142.9
1-N-Hydroxyethylflurazepam	150	133.3
Medazepam	150	133.3
Flurazepam	190	105.3
Lormetazepam	200	100.0
Clonazepam	210	95.2
Oxazepam	250	80.0
Clobazam	260	76.9
Clotiazepam	380	52.6
Lorazepam	600	33.3
Bromazepam	630	31.7
7-Aminoflunitrazepam	930	21.5
Demoxepam	1600	12.5
7-Aminoclonazepam	5300	3.8
Norchlordiazepoxide	4500	4.4
Chlordiazepoxide	3300	6.1
Clorazepate	*	
Clorazepate degrades rapidly in stomach acid to nordiazepam. Nordiazepam hydroxylates to oxazepam.		

Cocaine Metabolite

Compounds	Concentration equivalent to 300 ng/mL Cutoff	% Cross-Reactivity
Benzoyllecgonine	300	100
Ecgonine	7000	4.3
Cocaine	40000	0.8

Fentanyl

Norfentanyl (Major Metabolite)

Compounds	Concentration Tested (ng/mL)	% Cross-Reactivity
Norfentanyl (Major Metabolite)	15	7

Other Metabolites and Structural Analogs of Fentanyl

Compounds	Concentration equivalent to 1 ng/mL Cutoff	% Cross-Reactivity
Acetyl fentanyl	1.1	90.91
Isobutyryl fentanyl	1.1	90.91
ω -1-Hydroxyfentanyl	1.2	83.33
Acrylfentanyl	1.3	76.90
Butyryl fentanyl	1.4	71.43
Furanyl fentanyl	1.5	66.67
Para-fluoro fentanyl	1.5	66.67
Ocfentanil	1.6	62.50
4-Fluoro-isobutyryl fentanyl	1.9	52.63
Para-fluorobutyryl fentanyl (p-FBF)	1.9	52.63
Valeryl fentanyl	2.3	43.48
β -hydroxyfentanyl	9.5	10.53
Acetyl norfentanyl	12.1	8.26
(\pm) β -hydroxythiofentanyl	32.7	3.06
(\pm)-3-cis-methyl fentanyl	144.1	0.69
Carfentanil	448.2	0.22
Despropionyl fentanyl (4-ANPP)	471.8	0.21
Sufentanil	2,362	0.04
Remifentanil	10,000	<0.01
Norcarfentanil	38,196	0.003
Alfentanil	100,000	<0.001

Methadone

Compounds	Concentration equivalent to 300 ng/mL Cutoff	% Cross-Reactivity
Methadone	300	100

Opiates

Compounds	Concentration equivalent to 300 ng/mL Cutoff	% Cross-Reactivity
Codeine	102	294.1
Ethylmorphine	240	125.0
Hydrocodone	247	121.5
Dihydrocodeine	291	103.1
Morphine	300	100.0
6-Acetylmorphine	435	69.0
Levorphanol	480	62.5
Hydromorphone	498	60.2
Morphine-3-Glucuronide	626	47.9
Nalorphine	2130	14.1
Oxycodone	3340	9.0
Levallorphan	3740	8.0
Oxymorphone	9300	3.2
Meperidine	>15000	<2
Naloxone	360000	0.1

Oxycodone

Compounds	Concentration equivalent to 100 ng/mL Cutoff	% Cross-Reactivity
Oxycodone	100	100
Oxymorphone	97	103

Phencyclidine (PCP)

Compounds	Concentration equivalent to 25 ng/mL Cutoff	% Cross-Reactivity
Phencyclidine (PCP)	25	100
4-Phenyl-4-piperidinocyclohexanol	32	78.1
1-[1-(2-Thienyl)-cyclohexyl]piperidine (TCP)	37	67.6
1-(1-Phenylcyclohexyl)morpholine (PCM)	41	61.0
1-(1-Phenylcyclohexyl)pyrrolidine (PCPy)	54	46.3
1-[1-(2-Thienyl)-cyclohexyl]morpholine (TCM)	80	31.3
1-[1-(2-Thienyl)-cyclohexyl]pyrrolidine (TCPy)	83	30.1
N,N-Diethyl-1-phenylcyclohexylamine (PCDE)	234	10.7
1-(4-Hydroxypiperidino)phenylcyclohexane	420	6.0
Mesoridazine	50000	0.1
Meperidine	67000	<0.1
Dextrorphan	97000	<0.1
Dextromethorphan	120000	<0.1

THC/Cannabinoids

Compounds	Concentration equivalent to 50 ng/mL Cutoff	% Cross-Reactivity
11-nor- Δ -9-THC-9-COOH	50	100
8- β -11-Dihydroxy- Δ -9 -THC	58	86
8- β -Hydroxy- Δ -9 -THC	68	74
11-Hydroxy- Δ -8 -THC	67	75
11-Hydroxy- Δ -9 -THC	77	65
9-Carboxy-11-nor- Δ -9-THC-glucuronide	95	53