

Beckman AU Drug Screen Cutoffs and Calibrator Identification

Drug Class	Calibrator	Cutoff (ng/mL)
Amphetamines	d-Methamphetamine	1000
Barbiturates	Secobarbital	200
Benzodiazepines	Lormetazepam	200
Buprenorphine	Norbuprenorphine	5
Cocaine Metabolite	Benzoylcegonine	300
Fentanyl	Fentanyl	1.0
Methadone	Methadone	300
Opiates	Morphine	300
Oxycodone	Oxycodone	100
Phencyclidine (PCP)	Phencyclidine (PCP)	25
Tetrahydrocannabinol (THC)	11-nor- Δ -9-THC-9 Carboxylic Acid	50