

Abbott Alinity Drug Screen Cutoffs and Calibrator Identification

| Drug Class | Calibrator | Cutoff (ng/mL) |
|----------------------------|-------------------------------------------|-----------------------|
| Amphetamines | d-Methamphetamine | 1000 |
| Barbiturates | Secobarbital | 200 |
| Benzodiazepines | Oxazepam | 200 |
| Buprenorphine | Buprenorphine | 5 |
| Cocaine Metabolite | Benzoylcegonine | 300 |
| Fentanyl | Fentanyl | 1.0 |
| Methadone | Methadone | 300 |
| Opiates | Morphine | 300 |
| Oxycodone | Oxycodone | 100 |
| Phencyclidine (PCP) | Phencyclidine (PCP) | 25 |
| Tetrahydrocannabinol (THC) | 11-nor- Δ -9-THC-9 Carboxylic Acid | 50 |